remel

x/pect® Giardia Control Kit

INTENDED USE

Remel Xpect® Giardia Control Kit is a test kit containing characterized positive and negative controls intended for use in qualitative procedures as external quality control material for Xpect® Giardia or Xpect® Giardia/Cryptosporidium test kits. These controls are intended to monitor gross systematic failures only.

REAGENTS AND MATERIALS SUPPLIED

(approximately 15 tests)

- Giardia Positive Control (1 x 2 ml):**
Control material derived from 10% Formalin or SAF fixed stool samples of human origin, which are positive for *Giardia*
- Giardia/Cryptosporidium Negative Control (1 x 2 ml):**
Control material derived from 10% Formalin or SAF fixed stool samples of human or animal origin, which are negative for *Giardia* and *Cryptosporidium*
- Instructions For Use (IFU)**

MATERIALS REQUIRED BUT NOT SUPPLIED

- Xpect® Giardia Test Kit
(REF R2450020), or
- Xpect® Giardia/Cryptosporidium Test Kit
(REF R2450520)

STORAGE

Store all components at 2-8°C. Do not freeze or overheat. This product should not be used beyond the expiration date printed on the outer package label.

PRECAUTIONS

- For *In Vitro* Diagnostic Use.
- Refer to the IFU for Xpect® Giardia or Xpect® Giardia/Cryptosporidium test kit for specific precautions related to the test.
- Controls should be considered potentially hazardous and handled in the same manner as an infectious agent.

PROCEDURE

It is recommended that external positive and negative controls be run with each new test kit lot number or shipment of product. Each laboratory should follow their local, state, and/or federal regulations or accreditation requirements for quality control testing.

- Allow the controls to reach room temperature (20-25°C) prior to testing. The controls are ready to use. No dilution is required.
- Process controls in accordance with the Procedure section contained in the most current IFU for Xpect® Giardia or Xpect® Giardia/Cryptosporidium test kit.
- Use four (4) full drops of the corresponding control (positive or negative) in place of a stool specimen in the procedure. Four drops from the control dropper equals approximately 0.1 ml.
- If controls do not perform as expected, do not report patient results.

INTERPRETATION

The test should produce positive and negative results with the corresponding control as outlined in the Xpect® test kit IFU. The test should be repeated if invalid results are obtained.

PACKAGING

REF R2450030Xpect® Giardia Control Kit

Symbol Legend

REF	Catalog Number
IVD	In Vitro Diagnostic Medical Device
LAB	For Laboratory Use
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
LOT	Batch Code (Lot Number)
	Use By (Expiration Date)
EC REP	European Authorized Representative

Xpect® is a registered trademark of Thermo Fisher Scientific and its subsidiaries.

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